



MedLog Solutions, LLC



Clarity Platinum Series



The Clarity Platinum Urine Analyzer is the first urine analyzer to be FDA Approved since the 2018 guidelines were created. The analyzer has cutting edge technology to ensure compliance and accuracy. It is also HL7 Compliant and Made in the USA.

<p>Product Code: CLA-PLTUAPROMO5</p> <p>Contents:</p> <ul style="list-style-type: none"> ✓ Clarity Platinum Urine Analyzer ✓ 5 Bottles CLA-10P ✓ Barcode Reader ✓ 12 Month Warranty 	<p>Product Code: CLA-10P</p> <p>Contents:</p> <ul style="list-style-type: none"> ✓ Tests for: Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, and Glucose ✓ 100 Strips per bottle 	<p>Product Code: CLA-UHCRL25</p> <p>Contents:</p> <ul style="list-style-type: none"> ✓ Clarity Urinalysis Controls 1X25 mL -Negative 1 X 25 mL -Positive ✓ For use with CLA-10P & CLA-MAC 	<p>Product Code: DTG-UAPPR</p> <p>Contents:</p> <ul style="list-style-type: none"> ✓ Clarity Thermal Printer Paper, 4 rolls per pack <hr/> <p>Product Code: DTG-UASPPR</p> <p>Contents:</p> <ul style="list-style-type: none"> ✓ Clarity Thermal Sticker Printer Paper Rolls, 4 rolls per pack



MedLog Solutions, LLC
A Service-Disabled Veteran-Owned Small Business
Customer Service (844) 590-8787
Orders@MedLogSolutions.com

www.MedLogSolutions.com

Duns: 080705380
Cage Code: 7VWG2
UEI: XEKUDLKN3HZ4
FSS/GSA Contract:
36F79723D0174
Federal Tax ID:
82-1756285